


# MEDITECH Expanse TIP SHEET

## CRRT-P Cont Renal Replacement

The CRRT-P Cont Renal Replacement (CRRTPRISMA.C.OM) order was updated to account for a new fluid group response option (Phoxillum BK22K 4/0).

These changes were made in collaboration with the Corporate Service Lines and Provider Product teams



The image displays three screenshots of the CRRT-P Cont Renal Replacement order form, illustrating the changes made to the interface. Each screenshot shows a different field highlighted in red, indicating the previous options available when selecting "Type of Therapy" CVVH group response. The highlighted fields are:

- Pre blood pump (White scale)** (highlighted in the first screenshot)
- Post filter replacement fluid (Green scale)** (highlighted in the second screenshot)
- Pre filter replacement fluid (Purple scale)** (highlighted in the third screenshot)

These are the previous options that were available when selecting "Type of Therapy" CVVH group response. The following are highlighted

# MEDITECH Expanse TIP SHEET

## CRRT-P Cont Renal Replacement

	<p>These are the previous options that were available when selecting “Type of Therapy” CVVHD group response. The following are highlighted</p>
	<p>These are the previous options that were available when selecting “Type of Therapy” CVVHF group response. The following are highlighted</p>


# MEDITECH Expanse TIP SHEET

## CRRT-P Cont Renal Replacement

<p>CRRT-P Cont Renal Replacement</p> <ul style="list-style-type: none"> <li>* Blood flow rate</li> <li>* Type of Therapy: CVVHDF</li> <li>* Pre/Post filter replacement fluid (Purple scale)</li> <li>* Replacement solution rate (Purple scale)</li> <li>* Dialysate (Green scale)</li> <li>* Dialysate rate (Green scale)</li> <li>Pre blood pump (White scale)</li> </ul>	
<p>CRRT-P Cont Renal Replacement</p> <ul style="list-style-type: none"> <li>replacement fluid (Purple scale)</li> <li>* Replacement solution rate (Purple scale)</li> <li>* Dialysate (Green scale)</li> <li>* Dialysate rate (Green scale)</li> <li>Pre blood pump (White scale)</li> <li>Pre blood pump rate (White scale)</li> <li>* Net fluid removal (mL/hr)</li> <li>Anticoagulation: Yes No</li> </ul>	
<p>CRRT-P Cont Renal Replacement</p> <ul style="list-style-type: none"> <li>* Access type</li> <li>* Filter</li> <li>* Blood flow rate</li> <li>* Type of Therapy: CVVH</li> <li>* Pre filter replacement fluid (Purple scale)</li> <li>* Replacement solution rate (Purple scale)</li> <li>* Post filter replacement fluid (Green scale)</li> <li>Pre blood pump (White scale)</li> </ul>	<p>This is the new option that is available when selecting "Type of Therapy" CVVH group response. The following are highlighted</p>

# MEDITECH Expanse TIP SHEET

## CRRT-P Cont Renal Replacement



The image displays two screenshots of the 'CRRT-P Cont Renal Replacement' form. The top screenshot shows the 'Pre blood pump rate (White scale)' dropdown menu open, with 'Phoxillum B22K 4/0' highlighted by a red box. The bottom screenshot shows the 'Pre blood pump rate (White scale)' dropdown menu open, with 'Phoxillum B22K 4/0' highlighted by a red box. Both screenshots show the following fields:

- \* Type of Therapy: CVWH
- \* Pre filter replacement fluid (Purple scale): [Empty]
- \* Replacement solution rate (Purple scale): [Empty]
- \* Post filter replacement fluid (Green scale): [Empty]
- Pre blood pump (White scale): [Empty]
- Pre blood pump rate (White scale): [Empty]
- \* Net fluid removal (mL/hr): [Empty]


# MEDITECH Expanse TIP SHEET

## CRRT-P Cont Renal Replacement

<p>CRRT-P Cont Renal Replacement</p> <p>Blood flow rate</p> <p>Type of Therapy: CVVHD</p> <p><b>Dialysate (Green scale)</b>: [Dropdown menu open with 'Phoxillum B22K 4/0' highlighted in red]</p> <p>Dialysate rate (Green scale)</p> <p>Pre blood pump (White scale)</p> <p>Pre blood pump rate (White scale)</p> <p>Net fluid removal (mL/hr)</p> <p>Anticoagulation: Yes No</p> <p>Blood priming: Yes No</p>	<p>This is the new option that is available when selecting "Type of Therapy" CVVHD group response. The following are highlighted</p>
<p>CRRT-P Cont Renal Replacement</p> <p>Blood flow rate</p> <p>Type of Therapy: CVVHD</p> <p>Dialysate (Green scale)</p> <p>Dialysate rate (Green scale)</p> <p><b>Pre blood pump (White scale)</b>: [Dropdown menu open with 'Phoxillum B22K 4/0' highlighted in red]</p> <p>Pre blood pump rate (White scale)</p> <p>Net fluid removal (mL/hr)</p> <p>Anticoagulation: Yes No</p> <p>Blood priming: Yes No</p>	<p>This is the new option that is available when selecting "Type of Therapy" CVVHDF group response. The following are highlighted</p>
<p>CRRT-P Cont Renal Replacement</p> <p>Blood flow rate</p> <p>Type of Therapy: CVVHDF</p> <p><b>Pre/Post filter replacement fluid (Purple scale)</b>: [Dropdown menu open with 'Phoxillum B22K 4/0' highlighted in red]</p> <p>Replacement solution rate (Purple scale)</p> <p>Dialysate (Green scale)</p> <p>Dialysate rate (Green scale)</p> <p>Pre blood pump (White scale)</p> <p>Pre blood pump rate (White scale)</p>	<p>This is the new option that is available when selecting "Type of Therapy" CVVHDF group response. The following are highlighted</p>

# MEDITECH Expense TIP SHEET

## CRRT-P Cont Renal Replacement



The image displays two screenshots of the MEDITECH CRRT-P Cont Renal Replacement form. The top screenshot shows the 'Dialysate (Green scale)' dropdown menu open, with 'Phoxillum B22K 4/0' selected. The bottom screenshot shows the 'Pre blood pump (White scale)' dropdown menu open, also with 'Phoxillum B22K 4/0' selected. Both dropdown menus list other options: PrismaSOL BGK 4/2.5, PrismaSOL BK0/0/1.2, and Other.

**Top Screenshot:**

- CRRT-P Cont Renal Replacement
- Dialysate flow rate
- Type of Therapy: CVVHDF
- Pre/Post filter replacement fluid (Purple scale)
- Replacement solution rate (Purple scale)
- Dialysate (Green scale):** Phoxillum B22K 4/0
- Dialysate rate (Green scale)
- Pre blood pump (White scale): Phoxillum B22K 4/0
- Pre blood pump rate (White scale)

**Bottom Screenshot:**

- CRRT-P Cont Renal Replacement
- Replacement solution rate (Purple scale)
- Dialysate (Green scale)
- Dialysate rate (Green scale)
- Pre blood pump (White scale):** Phoxillum B22K 4/0
- Pre blood pump rate (White scale)
- Net fluid removal (mL/hr)
- Anticoagulation: Yes No