


Expanse TIP SHEET

Wound Image Capture




Expanse Mobile POC can add a wound image to a patient's chart.

After locating the patient in POC, use the More tab at the bottom to expand the menu and select the Note icon.

Select Add Note

Scroll down the list to find Wound Image Note.




Add the Location and Laterality information.

Select the Camera icon.
Features:

- Turn camera for full length image. System will automatically adjust the view.
- Ability to retake photos.

****Follow your facilities wound capture protocol for camera distance and measurements. It is best practice to take the image 6in from the wound.**



Select in the white space under the picture to add verbiage about the image.

Multiple images can be added here by selecting the camera icon again.

Select Save when finished.

Use the drop-down arrow to edit after saving if needed. There will be a record of the edit.

***It does take several minutes for the image to load in the patient's chart.

Software Tools

Home

Chart

Document

Library

Discharge

Diagnostics

History & Problems

Summary

Provider Notes

Administrative

Activity

Nurse/Allied Health

Other Clinical

Flowsheets

Nurse/Allied Health

All Activity | All Time | All Clinician Specialties | All Subcategories | Selected Visits

Nursing

Clinical Rehab Specialist

Hand Off


Occupational Therapy

Wound Care

DATE	DESCRIPTION	AUTHOR	TYPE
Today 16:37	Wound Care Note right leg Recorded: 06/19/24 16:36		Wound Care
Today 16:35	Wound Care Note left elbow Recorded: 06/19/24 16:37		Wound Care
06/13/24 10:47	Wound Care Note healing Recorded: 06/13/24 10:48		Wound Care
06/13/24 10:46	Wound Care Note		Wound Care

When viewing the patient's chart in Expanse Desktop, notice the Note icon containing the image.

***Always be sure to document all wound documentation within your wound assessment intervention(s).




After selecting the Note, it will open in an overlay.


To view "ALL WOUND CARE" or "ALL NOTES" for this patient use the tabs at the top.

For better viewing, right-click on the image to open in a new browser tab.

To view images side by side, pull the newly opened browser windows out by clicking it and dragging.



Use the three dots on the browser bar to access the Zoom button to zoom in on the image. This will affect all open browsers so be sure to readjust after viewing the image.






MEDITECH Expanse TIP SHEET

Febrile Neutropenia Protocol +

The Febrile Neutropenia Protocol + (12212810) order was updated to account for three new queries.

These changes were made in collaboration with the Corporate Service Lines and Provider Product teams

	This is the previous query for the Febrile Neutropenia Protocol + order
	Primary Antibiotic query was added as required – users to select from the group response options
	Additional Antibiotic query was added as Yes/No query – to determine if additional antibiotic is needed

MEDITECH Expanse TIP SHEET

Febrile Neutropenia Protocol +

Febrile Neutropenia Protocol +	
☒ .As Ordered	P New 05/07/26 11:54 ☆
* Provider	Hospitalist01,Provider
* Source	Intra-Procedure
* Primary Antibiotic to order for protocol use	
Additional Antibiotic to order for protocol use?	Yes No
* Additional Antibiotic to order for protocol use	Vancomycin Tobramycin
☒ .As Ordered P	☆

Upon selecting Yes to the 2nd query – the 3rd query will appear and will be required



If No is selected for the 2nd query it will remain suppressed

MEDITECH Expanse TIP SHEET

Measure CVP

The Measure CVP (MEASURECVP.C.OM) order was updated to add two new query options that are only available for patients 18 years or younger.

These changes were made in collaboration with the Corporate Service Lines and Provider Product teams

	For Measure CVP, this patient is 48 years old, and the two additional queries have been suppressed.
	For the updated order this patient is 13 years old and the query displays for the user to fill in.